

# Intermediate Fridge Range

COMBINING SUPERB  
PERFORMANCE AND QUALITY  
WITH GREAT VALUE.

## Key Features:

- QCPP Compliant
- Independently tested by the Pharmacy Guild Of NSW
- Self-closing glass door with inbuilt lock
- 24-month comprehensive warranty
- Rated to 38 Degrees



Medi Guard® 401 Fridge  
VMS Code: VACCMED4



Medi Guard® 601 Fridge  
VMS Code: VACCMED6

Medi Guard® 1001 Fridge  
VMS Code: VACCMED10



## Specifications:

Model	VACCMED4	VACCMED6	VACCMED10
<b>Dimensions External (W x D x H) mm</b>	606 x 620 x 1940	700 x 790 x 2058	11130 x 720 x 2023
<b>Gross Litres</b>	380L	600L	1000L
<b>Compliance</b>	QCPP and Cold Chain compliant + Pharmacy Guild Tested	TBC	QCPP and Cold Chain compliant + Pharmacy Guild Tested
<b>Warranty</b>	24 months parts and labour	24 months parts and labour	24 months parts and labour
<b>Weight</b>	75Kg	110Kg	142Kg
<b>Operating range</b>	2.C to 8.C	2.C to 8.C	2.C to 8.C
<b>Shelves</b>	7	6	12
<b>Lock</b>	Self closing lockable doors	Self closing lockable doors	Self closing lockable doors
<b>Alarm</b>	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm	High / Low and Door Open Audible and Visual Alarm
<b>Power</b>	10 AMP	10 AMP	10 AMP
<b>Monitoring</b>	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts	Medi Log 2 Data Logger Audible and Visible Alerts
<b>Cooling</b>	Fan Assisted	Fan Assisted	Fan Assisted
<b>Body Type</b>	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External	White Powder Coated Steel Internal and External
<b>Power Consumption</b>	4.935KWH (Over 24hrs)	6.209KWH (Over 24hrs)	9.903KWH (Over 24hrs)