



The ALLEVYN Difference:
The freedom to move

ALLEVYN Gentle Border Lite – a light and conformable¹ dressing that is ideal for skin tears and hard to dress² areas

 **smith&nephew**
ALLEVYN[®]
GENTLE BORDER
LITE
Wound Dressings

ALLEVYN[◇] Gentle Border Lite – the freedom to move

That's the ALLEVYN difference

Providing greater conformability*¹ making it the ideal dressing for skin tears. The ALLEVYN Gentle Border Lite range has a soft silicone gel adhesive for skin that is particularly fragile or sensitive.



A thin dressing that is kind to fragile skin

- Conformable¹ and light
- Contains a thinner foam core
- Easily molds to body contours and hard to dress areas²
- Easy to apply – can be lifted and repositioned without losing its adherent properties³
- Showerproof dressing²
- Minimal pain on removal³



Experience the ALLEVYN difference

Order Information		
Code	Size & Shape	Pkg
66800833	5cm x 5cm Square	10 per Box
66800836	5.5cm x 12cm Rectangle	10 per Box
66801041	8cm x 8.4cm Multisite	10 per Box
66801013	8.6cm x 7.7cm Oval	10 per Box

*In one direction

References

1. SMTL Report 10/3328/1, Conformability of wound dressings, August 2010
2. Smith & Nephew Data on File report OR-DOF/011
3. Smith & Nephew Data on File report OR-DOF/001

Smith & Nephew Pty Ltd
Australia
T +61 2 9857 3999
F +61 2 9857 3900

Smith & Nephew Ltd
New Zealand
T +64 9 820 2840
F +64 9 820 2841

[◇]Trademark of Smith & Nephew
All trademarks acknowledged
AL-0070-0111-anz V1 SN15000 REV0 (06/20)

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.